

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name: Taurine - ≥98%

Product Number: W381306 CAS Number: 107-35-7

MDL: MFCD00008197
Formula: C2H7NO3S
Formula Weight: 125.15 g/mol

TEST

Appearance (Color)

Appearance (Form)
Infrared spectrum

Purity (Titration by NaOH)

Total Decenoic Acids

Solubility (Turbidity)

c = 5%; Water

Solubility (Color)

Arsenic (As)
Cadmium (Cd)
Mercury (Hg)
Lead (Pb)

Expiration Date Period

5 Years

Specification: PRD.1.ZQ5.10000011355

Specification

White

Powder or Crystals Conforms to Structure

98.0 - 102.0 %

Clear

Colorless

3 ppm1 ppm1 ppm10 ppm

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.