


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
AquaMetProduct Number: **901759**
CAS Number: 1414707-08-6Formula: C₃₉H₅₅Cl₃N₄ORu
Formula Weight: 803.21 g/mol
Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	Green to Dark Green to Dark Beige
Appearance (Form)	Powder
¹ H NMR Spectrum	Conforms to Structure
Titration with AgNO ₃ After Oxygen Combustion	89.5 - 110.5 %
ICP Major Analysis Confirms Ruthenium Component	Confirmed

Specification: PRD.0.ZQ5.10000101817

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.