

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Chromotropic acid disodium salt dihydrate technical grade

Product Number:

126225

CAS Number:

5808-22-0

MDL:


MFCD00150612

Formula:

C₁₀H₆Na₂O₈S₂ · 2H₂O

Formula Weight:

400.29 g/mol



TEST

Specification

Appearance (Color)
Grey to Brown to Olive Green

Conforms to Requirements

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Purity (Titration by NaOH)

≥ 78.5 %

Water (by Karl Fischer)

Remarks:

Specification Date : 12/06/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.