

sigma-aldrich.com

3050 Spruce Street, Saint Louis, MO 63103, USA Website: www.sigmaaldrich.com Email USA: techserv@sial.com Outside USA: eurtechserv@sial.com

Product Specification

Product Name: Tin - foil, thickness 0.5 mm, 99.998% trace metals basis

| Product Number: |
|-----------------|
| CAS Number: |
| MDL: |
| Formula: |
| Formula Weight: |

265756 7440-31-5 MFCD00133862 Sn 118.71 g/mol

Sn

| TEST | Specification |
|--|--------------------------|
| | Conforme |
| Appearance (Color) Silver | Conforms |
| | Conforms |
| Appearance (Form) | Contornis |
| Foil | Contorne |
| Thickness | Conforms |
| 0.5 mm | |
| ICP Major Analysis | Confirmed |
| Confirms Tin Component | |
| Purity | Meets Requirements |
| 99.998% Based On Trace Metals Analysis | |
| Trace Metal Analysis | <pre>< 25.0 ppm</pre> |

Specification: PRD.3.ZQ5.10000021785

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.