


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2,3-Pyrazinedicarboxamide - 95%

Product Number: P56151
CAS Number: 6164-78-9
MDL: MFCD00006133
Formula: C₆H₆N₄O₂
Formula Weight: 166.14 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
White to Tan	
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Carbon	41.0 - 45.8 %
Nitrogen	31.8 - 35.6 %

Specification: PRD.1.ZQ5.10000036582

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.