

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

1,4-Cyclohexanedimethanol vinyl ether, mixture of cis and trans - 85%, technical grade

Product Number: 409952
CAS Number: 114651-37-5
MDL: MFCD00191966
Formula: C₁₀H₁₈O₂
Formula Weight: 170.25 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless to Yellow
Appearance (Form)	Liquid
Proton NMR spectrum	Conforms to Structure
GC Purity (sum of isomers)	≥ 83.0 %A

Specification: PRD.0.ZQ5.10000054376

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.