


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Phenacyltriphenylphosphonium bromide - technical grade, ≥90%

Product Number: 151335
CAS Number: 6048-29-9
MDL: MFCD00011920
Formula: C₂₆H₂₂BrOP
Formula Weight: 461.33 g/mol



TEST	Specification
Appearance (Color)	White to Yellow
Appearance (Form)	Powder or Solid or Crystals
Infrared spectrum	Conforms to Structure
Titration with AgNO ₃	90.0 - 110.0 %
After Oxygen Combustion	
Carbon	60.9 - 74.5 %

Specification: PRD.0.ZQ5.10000078082

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.