

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Sulfanilamide – Vetec™ reagent grade, 98%

Product Number:

V900101

CAS Number:


63-74-1

Formula:

C₆H₈N₂O₂S

Formula Weight:

172.20 g/mol

**TEST****Specification**

Appearance (Colour)

White to Off-White

Appearance (Form)

Powder

Solubility (Colour)

Colorless to Faint Yellow

Solubility (Turbidity)

Clear

50 mg/mL, 0.5M HCl

Purity (Titration by HCL04)

≥ 97 %

Specification: PRE.1.ZQ5.10000002590

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.