


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Cyanomethyl dodecyl trithiocarbonate - 98% (HPLC)

Product Number: **723029**
CAS Number: 796045-97-1



Formula: C₁₅H₂₇NS₃
Formula Weight: 317.58 g/mol
Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Color) White to Orange to Tan	Conforms to Requirements
Appearance (Form) Powder or Crystals	Powder or Solid
Proton NMR spectrum	Conforms to Structure
Carbon	55.3 - 58.1 %
Nitrogen	4.3 - 4.5 %
Purity (HPLC)	≥ 97.5 %

Specification: PRD.1.ZQ5.10000009658

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.