

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Tetrahydrofurfuryl methacrylate - contains 75 ppm HQ as inhibitor, 900 ppm MEHQ as inhibitor, 97%

Product Number:

**409456**

CAS Number:

2455-24-5

MDL:


MFCD00014515

Formula:

C<sub>9</sub>H<sub>14</sub>O<sub>3</sub>

Formula Weight:

170.21 g/mol

**TEST****Specification**

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 96.5 %

HQ as Inhibitor

50 - 100 ppm

MEHQ as Inhibitor

800 - 1000 ppm

Specification: PRD.1.ZQ5.10000020635

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.