

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti-Catalase antibody produced in mouse - ~2 mg/mL, clone CAT-505, purified immunoglobulin, buffered aqueous solution

Product Number: C0979

Storage Temperature: -20 °C

TEST Specification

Appearance (Form)
Appearance (Turbidity)
Appearance (Color)
IgG Concentration (OD 280nm)
WB-Cell Line/Tissue Extract
WB Working Concentration
Recommended Retest Period

4 Years

Specification: PRE.0.ZQ5.10000016487

Liquid Clear Colorless 2.0 - 2.5 mg/ml HepG2 Human Carcinoma 0.5 - 1.0 µg/ml

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.