

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Aluminum nitride - powder, 10 µm, ≥98%

Product Number:**241903**

CAS Number:

24304-00-5

MDL:

MFCD00003429

Formula:

AlN

Formula Weight:

40.99 g/mol

AlN**TEST****Specification**

Appearance (Color)

Grey

Appearance (Form)

Powder

X-Ray Diffraction

Conforms to Structure

Gravimetric Analysis

64.5 - 67.1 %

% Al

Average Particle Size

≤ 10.00 micron

Specification: PRD.0.ZQ5.10000097934

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.