

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

3-Chloropropylamine hydrochloride - 98%

Product Number:

142549

CAS Number:

6276-54-6

MDL:

MFCD00012913

Formula:

C₃H₈ClN · HCl

Formula Weight:

130.02 g/mol

**TEST****Specification**

Appearance (Color)

White to Tan

Appearance (Form)

Conforms to Requirements

Powder, Crystals, Crystalline Powder and/or Chunks

Infrared spectrum

Conforms to Structure

Titration with AgNO₃

97.5 - 102.5 %

After Oxygen Combustion

GC (area %)

≥ 97.5 %

Specification: PRD.2.ZQ5.10000004397

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.