

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

2,3-Dimethoxy-1,4-naphthoquinone - $\geq 99\%$, solid

Product Number:

D5439

CAS Number:

6956-96-3

MDL:

MFCD00184331

Formula:

C₁₂H₁₀O₄

Formula Weight:

218.21 g/mol

Storage Temperature:

-20 °C

TEST**Specification**

Appearance (Color)

Yellow

Appearance (Form)

Powder

Solubility (Color)

Yellow

Solubility (Turbidity)

Clear

10 mg/mL, DMSO

Proton NMR spectrum

Conforms to Structure

Purity (HPLC)

 $\geq 99\%$

Specification: PRD.1.ZQ5.10000039617

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.