


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Lithium carbonate - 99.997% trace metals basis

Product Number: 203629
CAS Number: 554-13-2
MDL: MFCD00011084
Formula: Li_2CO_3
Formula Weight: 73.89 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Titration with HCl % CO_3	79.2 - 83.3 %
ICP Major Analysis Confirms Lithium Component	Confirmed
Purity 99.997% Based On Trace Metals Analysis	Meets Requirements
Trace Metal Analysis	\leq 35.0 ppm

Specification: PRD.1.ZQ5.10000047333

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.