


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Benzylhydrazine dihydrochloride - 97%

Product Number: **B22852**
CAS Number: 20570-96-1
MDL: MFCD00012921
Formula: C₇H₁₀N₂ · 2HCl
Formula Weight: 195.09 g/mol



TEST	Specification
Appearance (Color) White to Yellow to Beige	Conforms to Requirements
Appearance (Form) Powder and/or Chunks	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Carbon	41.6 - 44.6 %
Nitrogen	13.9 - 14.9 %
Solubility (Turbidity) C = 50mg/ml, 95% Ethanol	Clear
Solubility (Color)	Colorless to Dark Yellow

Specification: PRD.1.ZQ5.10000012597

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.