


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Poly(4-vinylpyridine) – average M_w ~60,000

Product Number: 472344
CAS Number: 25232-41-1
MDL: MFCD00190588
Formula: (C₇H₇N)_n



TEST

Specification

Appearance (Color)
Appearance (Form)
Powder or Granular Powder
Infrared spectrum
Loss on Drying

White to Off-White
Conforms to Requirements
Conforms to Structure
≤ 5.0 %

Specification: PRD.0.ZQ5.10000031515

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.