

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Glass beads, unwashed - $\leq 106 \mu\text{m}$ (140-finer U.S. sieve)

Product Number:

G8893

TEST	Specification
Appearance (Color)	White to Light Grey
Appearance (Form)	Beads
Bead Size	Conforms
106 Microns and finer	

Specification: PRD.0.ZQ5.10000059072

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.