

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Ammonia solution 0.4 M in THF

Product Number:

718939

NH₃

TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Purity Titration by NaOH (M)	0.30- 0.50 M
Recommended Retest Period	-----
Six Months	

Remarks:

Specification Date : 09/22/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.