

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Muscarinic Acetylcholine Receptor (M3) antibody produced in rabbit - affinity isolated antibody, lyophilized powder

Product Number: **M0194**

MDL: MFCD06798471

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Formulation	Conforms
LYOPHILIZED FROM PHOSPHATE BUFFERED SALINE, pH 7.4, 1% BSA AND 0.05% SODIUM AZIDE	
Concentration	Conforms
0.6-0.8 MG/ML AFTER RECONSTITUTION	
Working Dilution	Conforms
By Direct Immunoblotting	
1:200 USING RAT BRAIN LYSATE	

Specification: PRD.0.ZQ5.10000043908

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.