


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Bromophenol Blue sodium salt

Product Number: B5525
CAS Number: 34725-61-6
MDL: MFCD00013793
Formula: C₁₉H₉Br₄NaO₅S
Formula Weight: 691.94 g/mol



TEST	Specification
Appearance (Color) Dark blue to black with a gray cast	Conforms
Appearance (Form)	Powder
Suitability Suitable for Use as a Tracking Dye in Nucleic Acid Gel Electrophoresis.	Suitable
Protease Detection	None Detected

Specification: PRD.1.ZQ5.1000009346

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.