


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Polyethylenimine, branched - average M_w ~270,000 by LS

Product Number: **904759**
CAS Number: 9002-98-6



Formula: H

TEST	Specification
Appearance (Color)	Colorless to Yellow
Appearance (Form)	Liquid
^1H NMR Spectrum	Conforms to Structure
Water (by Karl Fischer)	$\leq 12.0 \%$
Molecular Weight by LS	220000 - 320000

Specification: PRD.0.ZQ5.10000101251

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.