

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Poly(acrylic acid) partial sodium salt solution – average M_w ~2,000 by GPC, 61–65 wt. % in H₂O

Product Number:

192023


CAS Number:

9033-79-8

MDL:

MFCD00084394

Formula:

(C₃H₄O₂.C₃H₄O₂.Na)_x**TEST****Specification**

Appearance (Color)

Colorless to Pale Yellow

Appearance (Form)

Viscous Liquid

Infrared spectrum

Conforms to Structure

Brookfield Viscosity

400 - 1400 cps

pH

2.2 - 3.0

Total Solids

62.0 - 64.0 %

Specification: PRD.0.ZQ5.10000051323

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.