

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:


2-(Methylamino)propiofenone hydrochloride solution – drug standard, 1.0 mg/mL in methanol

Product Number:

M8536

MDL:

MFCD00012380



Storage Temperature:

2 - 8 °C

TEST

Specification

Concentration

0.95 - 1.05 mg/ml

Free Base

Infrared spectrum

Conforms to Structure

Expiration Dated Product

Specification: PRD.3.ZQ5.10000050013

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.