

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

MTT Formazan - powder

 Product Number:
 M2003

 CAS Number:
 57360-69-7

 MDL:
 MFCD00080731

 Formula:
 C18H17N5S

 Formula Weight:
 335.43 g/mol

() N S OH

TEST

Appearance (Form)
Carbon
Nitrogen
Purity (TLC)
Solubility (Color)

c = 1 mg/ml; Acetone
Dark Purple

Specification: PRD.1.ZQ5.10000013327

Specification

Powder 62.8 - 66.1 % 20.4 - 21.4 % > 97 %

Conforms to Requirements