

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:


Poly(ethyleneimine) solution average M_n ~1,200, average M_w ~1300 by LS, 50 wt. % in H₂O

Product Number:

482595

MDL:

MFCD00084427



TEST	Specification
Appearance (Color)	Colorless to Yellow
Appearance (Form)	Viscous Liquid
Infrared spectrum	Conforms to Structure
Concentration	48.0- 52.0 %
Viscosity sec (flow cup)	40- 100

Remarks:

Specification Date : 11/24/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.