

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Human IgG (y-chain specific), F(ab')2 fragment-FITC antibody produced in goat - affinity isolated antibody, buffered aqueous solution

Product Number: F1641

MDL: MFCD00162455

Storage Temperature: 2 - 8 ℃

TEST Specification

Appearance (Turbidity) Clear

Appearance (Color) Faint Yellow to Yellow

Appearance (Form) Liquid
F/P Molar Ratio 3.0 - 5.0

Prior to the addition of BSA

A280/A496 UV Absorption 1.0 - 1.5

Prior to the addition of BSA

Protein Content

Prior to the addition of BSA

Immunoflourescence (Direct) Minimum Dilution 1:32 using Human Peripheral Blood Lymphocytes

Specification: PRD.1.ZQ5.10000095987

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.