

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

N¹-(3-Trimethoxysilylpropyl)diethylenetriamine - technical grade

Product Number:

413348

CAS Number:

35141-30-1

MDL:


MFCD00069173

Formula:

C₁₀H₂₇N₃O₃Si

Formula Weight:

265.43 g/mol

**TEST****Specification**

Appearance (Color)

Colorless to Yellow

Appearance (Form)

Liquid

Proton NMR spectrum

Conforms to Structure

13C NMR Spectrum

Conforms to Structure

Titration by HCL

Typically > 85%

Specification: PRD.1.ZQ5.10000007618

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.