

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Silver – evaporation slug, diam. × L 0.6 cm × 1.2 cm, 99.99% trace metals basis

Product Number:

373249

CAS Number:

7440-22-4

MDL:

MFCD00003397

Formula:

Ag

Formula Weight:

107.87 g/mol

Ag

TEST	Specification
Appearance (Color) Silver	Conforms
Appearance (Form) Pellets or Slugs	Conforms
Size 0.6 cm Diameter x 1.2 cm Length	Conforms
ICP Major Analysis Confirms Silver Component	Confirmed
Trace Metal Analysis	≤ 150.0 ppm
Purity 99.99% Based On Trace Metals Analysis	Meets Requirements

Specification: PRD.2.ZQ5.10000005488

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.