


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Nickel(II) hexafluoroacetylacetonate hydrate - 98%

Product Number: **339709**
CAS Number: 207569-13-9
MDL: MFCD00150916
Formula: $C_{10}H_2F_{12}NiO_4 \cdot xH_2O$
Formula Weight: 472.79 g/mol



TEST	Specification
Appearance (Color) Green	Conforms to Requirements
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Carbon (anhydrous)	24.7 - 26.0 %
Titration with EDTA % Ni	≤ 12.2 %
Water (by Karl Fischer)	≥ 1.9 %

Specification: PRD.1.ZQ5.10000077770

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

