


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: S8876
CAS Number: 127-79-7
MDL: MFCD00023212
Formula: C₁₁H₁₂N₄O₂S
Formula Weight: 264.30 g/mol

**TEST****Specification**

| | |
|---------------------------------|---------------------------|
| Appearance (Color) | White to Off-White |
| Appearance (Form) | Crystalline Powder |
| Solubility (Color) | Colorless to Faint Yellow |
| Solubility (Turbidity) | Clear |
| 50 mg/mL, 1M NH ₄ OH | |
| Titration by NaNO ₂ | ≥ 99.0 % |
| Recommended Retest Period | ----- |
| 4 Years | |

Specification: PRD.0.ZQ5.10000048714

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.