

3050 Spruce Street, Saint Louis, MO 63103, USA  
Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)  
Email USA: [techserv@sial.com](mailto:techserv@sial.com)  
Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

**Product Specification**

Product Name:  
Monoclonal Anti-Collagen, Type IV antibody produced in mouse – clone COL-94, ascites fluid

Product Number: C1926

Storage Temperature: -20 °C

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Very Faint Yellow to Very Dark Yellow to Very Dark Brown
Protein by Biuret	15 - 55 mg/ml
Mouse IgG1 by RID	4.0 - 8.0 mg/ml
IHC - Tissue	Human Placenta or Human Tongue
IHC - Fixation	FFPE-Formalin Fixed Paraffin Embedded
IHC - Titer	1:500-1:1000
Recommended Retest Period	-----
5 Years	

Specification: PRE.0.ZQ5.10000015969

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.