

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Leukotriene B₄ - ~100 µg/mL in ethanol, ≥97%

Product Number:

L0517

CAS Number:

71160-24-2

MDL:

MFCD00135630

Formula:


C₂₀H₃₂O₄

Formula Weight:

336.47 g/mol

Storage Temperature:

-20 °C



TEST	Specification
Appearance (Color)	Colorless
Appearance (Turbidity)	Clear
Appearance (Form)	Liquid
µg/ml in EtOH	90 - 121
Purity (HPLC)	≥ 97 %

Specification: PRD.1.ZQ5.10000039066

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

