


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Gly-Arg-Gly-Asp-Thr-Pro - $\geq 97\%$ (HPLC)

Product Number: G5646
CAS Number: 108682-58-2
MDL: MFCD00076465
Formula: C₂₃H₃₉N₉O₁₀
Formula Weight: 601.61 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
1 mg/ml, H ₂ O	
Amino Acid Analysis	$\geq 60\%$
% Peptide Content	
Purity (HPLC)	$\geq 97\%$

Specification: PRD.0.ZQ5.10000045973

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.