


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-(2-Hydroxy-5-methylphenyl)benzotriazole - 97%

Product Number: 533203
CAS Number: 2440-22-4
MDL: MFCD00022903
Formula: C₁₃H₁₁N₃O
Formula Weight: 225.25 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Powder, Crystals or Crystalline Powder
Infrared spectrum
Purity (HPLC)

White to Yellow
Conforms to Requirements

Conforms to Structure
≥ 96.5 %

Specification: PRD.0.ZQ5.10000016623

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.