

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

L-Glutamic acid-1-¹³C - 99 atom % ¹³C, 98% (CP)

Product Number:

604968

MDL:


MFCD00144588

Formula:

C₅H₉NO₄

Formula Weight:

148.12 g/mol



TEST	Specification
Purity (HPLC)	≥ 98.0 %
Chiral Purity	≤ 1.0 %
¹³ C NMR Spectrum	Conforms
Proton NMR Spectrum	Conforms
¹³ C Enrichment	≥ 99.0
Appearance (Color)	White
Appearance (Form)	Powder

Specification: PRD.0.ZQ5.10000103468

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

