

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:


Ethylene glycol monopropyl ether – Vetec™ reagent grade

Product Number:

V900154

CAS Number:

2807-30-9



Formula:

C₅H₁₂O₂

Formula Weight:

104.15 g/mol

TEST	Specification
------	---------------

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Water (by Karl Fischer)

≤ 0.1 %

Specific Gravity

0.912 - 0.914

At 20 Degrees Celsius

Color Test

≤ 10 APHA

Purity

≥ 99.0 %

Specification: PRE.1.ZQ5.10000002618

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.