


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Dimethyloctadecyl[3-(trimethoxysilyl)propyl]ammonium chloride - 72%

Product Number: 435708
CAS Number: 27668-52-6
MDL: MFCD00054224
Formula: C₂₆H₅₈ClNO₃Si
Formula Weight: 496.28 g/mol



TEST	Specification
Appearance (Clarity) Assay Performed After Heating to 50 Degrees Celcius. May be Viscous.	Clear to Hazy
Appearance (Color)	Colorless to Yellow
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Supplier - Purity	69.1 - 74.9 %
Aegis Environmental Management Data	
Product of Supplier	Confirmed
Aegis Environmental Management	

Specification Date : 03/10/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.