


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Doxazosin mesylate - ≥97% (HPLC), powder

Product Number: D9815
CAS Number: 77883-43-3
MDL: MFCD00216023
Formula: C₂₃H₂₅N₅O₅ · CH₃SO₃H
Formula Weight: 547.58 g/mol



| TEST | Specification |
|---------------------|-----------------------|
| Appearance (Color) | White to Off White |
| Appearance (Form) | Powder |
| Proton NMR Spectrum | Conforms to Structure |
| Purity (HPLC) | ≥ 97 % |

Specification: PRD.1.ZQ5.10000068704

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.