

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Ficoll® 400 - Type 400-DL, lyophilized powder

Product Number: **F9378**  
CAS Number: 26873-85-8  
MDL: MFCD00081599

Formula Weight: 400 kg/mol

TEST	Specification
------	---------------

Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Very Faint Yellow
Solubility (Turbidity)	Clear to Slightly Hazy
100 mg/ml H <sub>2</sub> O	
µmhos/cm Conductivity	≤ 200
10% solution	
Recommended Retest Period	-----
5 years	
Note	-----
Ficoll is a registered trademark of Cytiva	

Specification: PRD.1.ZQ5.10000029764

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

