

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

β -Glucuronidase from Escherichia coli – $\geq 10,000,000$ units/g protein (30 min assay), recombinant, expressed in E. coli overproducing strain, lyophilized powder

Product Number:

G8295

CAS Number:

9001-45-0

MDL:

MFCD00130629

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Protein by Biuret

40 - 60 %

units/g protein

 ≥ 10000000

One Sigma or modified Fishman unit will liberate 1.0 microgram of phenolphthalein from phenolphthalein glucuronide per hour at 37 deg C at the pH 6.8 (30 min assay).

Recommended Retest Period

2 years

Specification: PRD.0.ZQ5.10000023006

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.