

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: 336807
CAS Number: 10025-83-9
MDL: MFCD00011063
Formula: Cl₃Ir
Formula Weight: 298.58 g/mol



TEST	Specification
Appearance (Color) Grey to Dark Brown or Dark Green	Conforms to Requirements
Appearance (Form) Crystals	Conforms to Requirements
Titration by AgNO ₃ After Oxygen Combustion	35 - 37 %

Specification: PRD.0.ZQ5.10000078087

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.