


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Benzyltriphenylphosphonium chloride - 99%

Product Number: B32807
CAS Number: 1100-88-5
MDL: MFCD00011913
Formula: C₂₅H₂₂ClP
Formula Weight: 388.87 g/mol



TEST

Specification

Appearance (Color)	White
Appearance (Form)	Powder or Crystals
Infrared Spectrum	Conforms to Structure
Titration with AgNO ₃	98.5 - 101.5 %
Solubility (Turbidity)	Clear
c = 25 mg/ml; CHCL ₃	
Solubility (Color)	Colorless
Purity (TLC)	≥ 99 %

Specification: PRD.1.ZQ5.10000017266

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.