

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

## **Product Specification**

Product Name:

Fructosyl-Amino Acid Oxidase from Corynebacterium sp. - recombinant, expressed in E. coli, lyophilized

powder, ≥0.45 units/mg protein

**Product Number:** F7296 CAS Number: 9037-65-4

MDL: MFCD03456889

Storage Temperature: -20 ℃

TEST Specification

units/mg solid > 0.45

One unit will produce 1.0 micromole of hydrogen peroxide per minute at pH 8.0 at 37 deg C.

Specification: PRD.0.ZQ5.10000074806

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.