


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: H5752
CAS Number: 68-96-2
MDL: MFCD00003659
Formula: C₂₁H₃₀O₃
Formula Weight: 330.46 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
EmM	16.0 - 16.8
EtOH	
Lambda Max	240 - 241 nm
Solubility (Color)	Colorless to Light Yellow
Solubility (Turbidity)	Clear
50 mg/mL, CHCl ₃	
Purity (HPLC)	≥ 95 %
Recommended Retest Period	-----
5 years	

Specification: PRD.0.ZQ5.10000038601

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.