

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Cadmium chloride - 99.99% trace metals basis

Product Number: 202908
CAS Number: 10108-64-2
MDL: MFCD00010916
Formula: CdCl₂
Formula Weight: 183.32 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Conforms to Requirements
Powder, Crystals or Solid and/or Chunks	
Compleximetric EDTA	60.1 - 62.5 %
% Cadmium	
ICP Major Analysis	Confirmed
Confirms Cadmium Component	
Purity	Meets Requirements
99.99% Based on Trace Metals Analysis	
Trace Metal Analysis	≤ 150.0 ppm

Specification: PRD.0.ZQ5.10000039723

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.