

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

3,3',5,5'-Tetramethylbenzidine (TMB) Liquid Substrate System for ELISA – peroxidase substrate

Product Number: **T0440**

Storage Temperature: 2 - 8 °C

TEST	Specification
pH	3.35 - 3.75
UV/VIS Absorbance at 285 nm	0.350 - 0.550
UV/VIS Absorbance at 650 nm	≤ 0.020
Expiration Dated Product	-----

Specification: PRD.4.ZQ5.10000012890

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.