

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Quinoxaline-6-boronic acid pinacol ester - 97%

Product Number:

708631

CAS Number:


1167418-13-4

Formula:

C₁₄H₁₇BN₂O₂

Formula Weight:

256.11 g/mol

**TEST****Specification**

Appearance (Color)

Colorless to Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Carbon

63.36 - 67.95 %

Nitrogen

10.55 - 11.32 %

Purity (HPLC)

≥ 96.5 %

Specification: PRD.0.ZQ5.10000045149

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.