

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

1',3'-Dihydro-8-methoxy-1',3',3'-trimethyl-6-nitrospiro[2H-1-benzopyran-2,2'-(2H)-indole] - 97%

Product Number:

433926

CAS Number:

1498-89-1

MDL:


MFCD00216608

Formula:

C<sub>20</sub>H<sub>20</sub>N<sub>2</sub>O<sub>4</sub>

Formula Weight:

352.38 g/mol



### TEST

### Specification

Appearance (Color)

Light Green to Very Dark Green and Brown  
to Very Dark Brown and Green-Brown and Orange Brown  
Powder

Appearance (Form)

Proton NMR Spectrum

Conforms to Structure

Purity (Titration by HCL04)

96.5 - 103.5 %

### Remarks:

Reviewing Chemist :

Specification: PRD.0.ZQ5.10000026585

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.