


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4'-Morpholinoacetophenone - 99%

Product Number: 119865
CAS Number: 39910-98-0
MDL: MFCD00006168
Formula: C₁₂H₁₅NO₂
Formula Weight: 205.25 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Faint Yellow to Yellow and Faint Brown to Light Brown
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Carbon	69.2 - 71.3 %
Nitrogen	6.5 - 7.1 %
Purity (GC)	≥ 98.5 %

Specification: PRD.0.ZQ5.10000040283

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.